

Twice as easy

Ambulatory actigraphy devices that collect and download continuous, objective long-term sleep/wake data



Actigraphy made easy

Our Actiwatch 2 and Actiwatch Spectrum ambulatory devices help make actigraphy easier to use, the devices more comfortable to wear, and the data acquisition and analysis functions less complex and time consuming.

The Actiwatch 2 incorporates a state-of-the-art ergonomic design, interchangeable watch bands, an ambient light sensor, and communication and software upgrades such as an automated rest algorithm. The Actiwatch Spectrum includes these features and a real watch face, off-wrist detection, color light recording and other advanced features. Now it's more efficient than ever for you to set up the Actiwatch system, gain user compliance and reliably collect and download sleep/wake data. That means it's now easier than ever to monitor a variety of patients who have a variety of sleep complaints. The Actiwatch 2 and Actiwatch Spectrum actigraphy systems combine improved aesthetics with advanced functionality to meet your needs.

Both devices incorporate high-quality biocompatible materials and well-engineered sensors (accelerometer, light, off wrist) to provide 24/7 wearability and accurate data collection. They have been thoroughly tested to be resistant to environmental conditions such as cold, heat, dust, perspiration and water. They also have been tested for drop resistance. Because of the unique strap material, both devices provide "wear it and forget it" comfort. Rugged durability helps enable the collection of complete data without any lapses.

The Actiwatch 2 and Actiwatch Spectrum can be used to analyze circadian rhythms, automically collect and score data for sleep parameters, and assess activity in any instance where quantifiable analysis of physical motion is desired.

Five modes

Both Actiwatch devices have five operational modes to help ensure battery efficiency and extend battery life.

- 1. Sleep
- 2. Standby
- 3. Data Collection
- 4. Post Collection
- 5. Communication




Actiware software

Actiware software helps make analyzing actigraphy data even easier than before.

Features include*:

- · Automatic rest algorithm
- Minor rest interval detection
- Charts up to seven channels of data on one actogram
- Supports legacy Actiwatch devices and ActiReaders
- Simultaneous data downloading from multiple devices
- Easy to set user defaults for common functions
- Automatic tabulation of sleep statistics for user- defined intervals
- Easy-to-use prompts for installation and communication to the Actiwatch devices
- Detailed how-to-use software tutorial website: www.learnactiware.com
- 21 CFR Part 11 compliant

*Version 5.5 or higher.




Wake/activity intervals are shown in black, sleep intervals in aqua, off wrist intervals in dark blue, event marks in dark blue triangles and light data in red, green, blue and yellow.



Actiwatch 2 special fetaures:

- Small, sleek and comfortable design to help enhance patient comfort and increase patient compliance
- Large memory (1 MBit non-volatile)
- Three data channels: photopic white light illuminance, activity, and event marking
- Ultrasonic welded case designed to provide moisture protection and increased resistance to environmental challenges (moisture, shock, dust, etc.)
- Faster communications rate
- Rechargeable battery via USB or wall current with 30-day life
- · Can be used with a hospital band
- Light-sensing records information about the amount and duration of ambient white light illuminance in units of lux
- Extra bands included with system




Wake/activity intervals are shown in black, sleep intervals in aqua, event marks in dark blue triangles and light data in yellow.



Actiwatch Spectrum special fetaures:

- Comfortable, aesthetic design to help enhance patient comfort and increase patient compliance
- · Off-wrist detector with real-time feedback on display
- Time and date indicators (like an actual wristwatch)
- Device status indicator (e.g., logging, not logging, on wrist, off wrist, and successful marker push)
- Seven data channels: photopic white light illuminance, red, green and blue light irradiance, activity, off wrist and event marking
- Moisture protection and increased resistance to environmental challenges (moisture, shock, dust, etc.)
- · Backlit display
- Infrared communication to the dock
- One-year battery life changeable at factory
- Large memory (2 MBit non-volatile)
- Faster communications and data-download rate
- · White light illuminance in units of lux
- Color light recording irradiance, in units of microwatts/cm², and photon flux, in units of photons/cm²/sec, in three color bands of the visual spectrum: red, green and blue
- Extra bands included with system



The free-living white light illuminance measurement performance as compared to a photometer is depicted, along with varying intensities of colored light irradiance measurements across common sources of light.

Philips Healthcare is part of Royal Philips Electronics

How to reach us

www.philips.com/healthcare Fax: +31 40 27 64 887

Asia

+852 2821 5888

Europe, Middle East, Africa +31 40 27 87246

Latin America +55 11 2125 0764

North America +1 425 487 7000 800 285 5585 (toll free, US only) Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668

Customer Service +1 724 387 4000 800 345 6443 (toll free)

Respironics Europe, Middle East, Africa +33 1 47 52 30 00

Respironics Asia Pacific +852 3194 2280

www.philips.com/respironics

For more information on the Actiwatch family of products, visit: actiwatch.respironics.com

Respironics, Actiwatch, ActiReader and Actiware are trademarks of Respironics, Inc. and its affiliates. All rights reserved. Windows is a registered trademark of Microsoft Corporation.



©2009 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

PrintTech SB 03/18/09 MCI 4102249 PN 1054519